
Panasonic®

Pore Cleanser

Operating Instructions

Model No. EH2591

Before operating this unit, please read these instructions completely.
Thank you for purchasing the Panasonic Pore Cleanser Model No. EH2591.

INDEX

ENGLISH Page 2

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Safety Instructions

- This personal care product is designed to remove excess oil and dirt on the face and nose area by application of a slight amount of suction. It is important that you exercise care when using this product in order to achieve satisfactory results and to reduce the risk of breaking capillaries beneath the skin surface. Depending on your complexion, you may notice purple spots or some skin discoloration as you first learn to use this product which will usually disappear within two to three days. This is dependent upon the length of time you apply the suction on a particular area of your face. If you are concerned about how your skin will react, we suggest that you test the product on a part of your face or neck that is not normally exposed. If discoloration occurs, discontinue use of the product.

WARNING

1. Do not hold either cup in one area, apply and quickly remove the cap.
2. Take care to limit the length of time you use the Pore Cleanser (up to one minute around the nose and about two minutes on the forehead and chin)

CAUTION

Purple spots or discoloration may appear on your skin.

- When using the nose cup, take care to use it only on the nose area. Do

not stop at any one point when using the cup.

- When using the forehead/chin cup, take care to use it only on the forehead and chin areas. Do not stop at any one point when using the cup.


Do not use on areas marked .

- Do not use the Pore Cleanser around your eyes, head, hairline border or on areas of thin skin.
- Do not use the Pore Cleanser on pimples, blemishes, sunburned or broken skin or any skin inflammations.
- If you have sensitive skin, take care not to use the Pore Cleanser on the same point more than twice.

If you use the Pore Cleanser soon after washing it with water, the suction force may be too strong, causing purple spots or scars to appear on your skin.



- After washing the Pore Cleanser, turn it on and leave it running for about 30 seconds. Then wipe the water from inside the Pore Cleanser before using again.

Operate for about 30 seconds.

Operating Instructions


Parts Identification

Installing the Battery


- * Make sure that your hands and the unit are dry during installation. Be careful not to allow water to get inside of the unit.

1. To remove the battery cover, insert the edge of a coin into the groove and rotate.

2. Insert one "AA" alkaline battery into the unit by aligning the \oplus and \ominus marks as shown in the illustration.

3. Reattach the battery cover and secure it by rotating the coin.

- Only use an "AA" alkaline battery. Other batteries may not provide sufficient power or perform for the same duration.
- Remove the batteries from the unit if it will not be used for a prolonged period of time.

Service life of the batteries

The typical service life for one "AA" alkaline battery is two months (based on five-minute periods of use, twice per week).


How to Use the Cups

Before using the Pore Cleanser for the first time, practice pulling the cup off the skin on the area of your upper arm.


Remove your finger from the pressure-release hole when pulling the unit away from your skin.

Do not pull away vertically from your skin.

Do not allow the unit to slide when pulling it away from your skin.


Treating Your Skin Using Suction

- For best results the pores of the skin should be well opened. For this reason it is a good idea to perform your skin care routine during or just after a bath or shower. The heat of the hot water will tend to open your pores.
 - We recommend that skin care with suction be performed two or three times a week.
 - When applying suction to your skin, place your finger over the pressure-release hole. Remove your finger from the pressure-release hole when pulling the unit away from your skin.
1. Remove any make-up from your face. Apply small amount of water on your face. Do not use when your skin is coated with cleansing lotion or facial soap.

2. Attach the cup socket, making sure that the nose cup is facing outward. Then turn the switch on.
3. Cover the pressure-release hole with your finger and move the Pore Cleanser slowly around the tip of the nose, as shown by the arrows. (See the picture on the next page) Remove your finger from the pressure-release hole when pulling the unit away from your skin.



Around the Tip of the Nose Using the nose cup

- ① Wet your skin slightly with water. Do not use cosmetic products or toiletries.
- ② Apply the Pore Cleanser around the tip of your nose by sliding it in the direction illustrated above.
- ③ Continue to treat your entire lower area of your nose for approx. one minute.

4. Fit the cup with the forehead/chin cup side facing outward, then turn the switch on again.
5. Covering the pressure-release hole down with your finger, move the Pore Cleanser over face from the center to the outside as shown by the arrows. (See the picture on the next page)

Forehead/chin cup

Forehead and Chin

- ① Wet your skin slightly with water. Do not use cosmetic products or toiletries.
 - ② Apply the Pore Cleanser around your forehead and chin by sliding it in the direction illustrated above.
 - ③ Continue to treat your entire forehead and chin for approx. two minutes.
- * Do not use the product on your cheeks, as the skin on this part of your face does not contain much sebum.

6. Rinse your face, and apply moisturizing lotion or any other skin conditioner that you normally use.

User-Maintenance Instructions

After using the Pore Cleanser both the main body and the cup can be washed using water.

1. Put some warm water in a washbowl. Do not use hot water (over 150°F (70°C))

Note: The pressure-release hole must be covered while cleaning.

* Do not use detergent.


2. Remove the cup. Turn the unit on, then soak the suction point in water, and clean for approximately 15 seconds while keeping the pressure-release hole covered.

3. Pull the main body out of the water, keeping the switch on with the pressure-release hole covered by a finger for about 30 seconds to allow water to drain from the main body.


- 4.** Wipe the main body and the cup with facial tissues or cotton swabs, taking particular care to clean carefully inside the cup.

How to Clean the Filter

Turn the switch off and remove any accumulated dirt using a toothbrush or similar tool like cotton swab.

Caution

- Not for use by children. Close supervision is necessary when this appliance is used near children or invalid.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.

- Never use while sleeping.
- Never insert an object into any opening.
- Do not use this appliance with a damaged or broken cup, as facial injury may occur.
- Do not submerge the hand unit for a long time.
Wash the hand unit in tap water. Do not use salt water, detergent or boiling water.
- The following people should consult a physician prior to using this product:
 - * People with acute illnesses, tumors or pregnant women
 - * People with skin types that are prone to allergic reactions
 - * People who have dermatitis, sunburnt or unusual skin conditions
 - * People with flesh wounds, rashes or abscesses

Complications could result from the use of this product by individuals who have the conditions described above.

- Never attempt to remodel, dismantle or repair the unit.
 - * Injuries could result from operating this unit under conditions for which it was not originally designed.

Troubleshooting Guide

Description	Probable cause	Action
Unit does not operate when turned on.	Battery is not installed.	Install battery.
	Battery cap is not properly secured.	Properly attach the battery cap.
No suction	The pressure hole is not being pressed.	Press the pressure-release hole.
Suction is weak.	Battery is weak.	Battery power is depleted; replace battery.
	The cup is not adhering to the skin.	Press cup against skin completely.
	The cup has not been attached to the unit properly.	Attach the cup properly.
Strong odor is being emitted.	The inside of the unit is not clean.	Wash out the unit; see "User-Maintenance Instructions".

Refrain from using the unit if you are unable to make it operate normally after reviewing the handling tips provided above.

Specifications

Power source : DC 1.5 V, one "AA" alkaline battery (not included)

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